

REMIDIO ANTERIOR IMAGING MODULE

AIM, or **Anterior Imaging Module** can digitize your slit lamp, enabling photo documentation and patient education.



new!
compatible
with all
iPhones from
6S onwards



bulk upload
with unlimited
cloud storage



redesigned
iPhone
adapter

high quality
beam splitter
optics

product specifications

model	AIM
smartphone	all iPhones from 6S onwards
camera resolution	as per iPhone
focal adjustment	manual
device weight without phone	900 grams
dimensions (w x l x h)	107 x 160 x 193mm



visit remidio.com

mail contactus@remidio.com

call +91 767 675 9999

