

CFA

C3 FIELD ANALYZER

Visual field screenings in less than 6 minutes with advanced algorithms



Perimetry on the go

Engineered by



Developed by





maintenance and
infrastructure free

product specifications

model name	C3FA
stimulus duration	200ms
stimulus interval	1300ms
reliability parameters	fixation loss, false positive
stimulus color	white on white perimetry
method of presentation	lcd screen
stimulus size	goldman III
general testing features	fixation monitoring via heijl-krakau method
screening test strategies	cfa supra threshold
threshold test strategies	cfa fast threshold strategy
screening test modes	multiple intensity
central patterns	10-2 24-2
fixation monitoring	heijl/krakau blind spot monitoring method
dimension	110 x 150 x 100
weight	~600 gms
max. temporal range	27 degrees
database management	cloud connectivity with cfa

visit remidio.com

mail contactus@remidio.com
info@alfaleustech.com

call +91 767 675 9999